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Bluetooth標準規格Ver. 5.0 (Bluetooth Low Energy)

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Startup Guide (this manual)

The Startup Guide describes the preparation procedure until injecting the odorant into the custom cartridges that can be purchased separately.

Help Guide (Web manual)

The Help Guide is an online manual that can be viewed via the Internet using a computer or smartphone.

For information about the Next-Gen Olfactometry System and how to use the custom cartridges (sold separately), refer to the following website.

https://rd1.sony.net/help/nos/nosdx1000/h_zz/

Safety Information

Manufacturer: This product has been manufactured by or on behalf of Sony Corporation, 1-7-1 Konan Minato-ku Tokyo, 108-0075 Japan.

EU and UK Importer/manufacturer's authorized representative is available at the following internet address:
<https://www.sony.net/nos-support-global/>

Communication specification

Frequency band: 2402 MHz – 2480 MHz

Max. EIRP (dBm): 6 dBm

Hereby, Sony Corporation declares that this equipment is in compliance with Directive 2014/53/EU and UK relevant statutory requirements.

The full text of the declaration of conformity is available at the following internet address:

<https://www.sony.net/nos-eu-compliance/>

**Disposal of old electrical & electronic equipment
(applicable in the European Union and other
countries with separate collection systems)**

This symbol on the product or on its packaging indicates that this product shall not be treated as household waste. Instead it shall be handed over to the appropriate collection point for the recycling of electrical and electronic equipment. By ensuring that this product is disposed of correctly, you will help to prevent potential negative consequences for the environment and human health, which could be caused by inappropriate waste handling of this product. The recycling of materials will help to conserve natural resources. For more detailed information about recycling of this product, please contact your local Civic Office, your household waste disposal service or the shop where you purchased the product.

Main unit

■ Caution

- Do not disassemble or modify the unit.
- Do not allow liquids or foreign objects to get inside the unit. If liquid should get inside, turn off the power immediately, unplug the power cord, and contact the dealer where the unit was purchased.
- Do not subject the unit to a strong impact. Doing so might cause an explosion, overheating or malfunction.
- Do not use anything other than dedicated custom cartridges (sold separately).
- Use the supplied AC adapter.
- Be careful not to pinch your fingers when replacing custom cartridges (sold separately).
- For hygienic reasons, replace the nose guide (page 5) with a new one for each user.

■ On placement

- Do not use the unit near water sources.
- Do not use, store or leave the unit at extremely high temperatures or in a place with high humidity.
- Place the unit on a flat and stable surface.

■ On AC adapter/cord

- Do not use the unit with a voltage exceeding the rating of the wall outlet. Exceeding the rating due to "octopus wiring" (overloading the outlet with too many plugs), etc. may cause fire due to heat generation.
- Insert the power plug firmly all the way to the base of the plug. Do not connect to a wall outlet that is loose even when plugged in all the way to the base.
- When unplugging the power cord, unplug it from the wall outlet first. Be sure to unplug it by holding the plug, not the cord.

Custom cartridge

■ Caution

- Do not modify the custom cartridge.
- Do not allow liquids or foreign objects other than odorants to get inside. If an abnormality occurs, immediately turn off the unit, unplug the power cord, and contact the dealer where the unit was purchased.
- Do not use anything other than dedicated custom cartridges (sold separately) and sub-cartridges (sold separately).
- Do not reuse a sub-cartridge (sold separately) that was injected with an odorant one time. Doing so may cause the odorant to leak.
- Do not open a sealed sub-cartridge (sold separately) that was injected with an odorant. Doing so may cause the odorant to leak.
- Follow the procedures in this manual and the Help Guide. Incorrect procedures may cause injury or damage.
- Perform the odorant injection workflow on a flat, stable surface.
- Be careful not to pinch your fingers when installing or removing a sub-cartridge (sold separately) from a custom cartridge (sold separately). Doing so may result in injury.
- Use the cartridge clasper (sold separately) to open or close the sub-cartridge.
- To use the torque wrench properly, read the operating instructions supplied with the torque wrench carefully.

■ On user's health

- If allergy symptoms appear, stop using the unit immediately.
- Examiners should stop using the unit immediately if the user feel unwell.
- If the user is pregnant or nursing, the examiner should use the unit with special attention.

- If dust or debris accumulates between the power plug and the wall outlet and moisture is absorbed, it may cause insulation degradation, resulting in a fire. Periodically unplug the power plug from the wall outlet to remove dust and debris.
- Unplug the power cord from the unit when moving the unit.
- For safety reasons, unplug the power cord from the wall outlet when the unit will not be used for an extended period of time.
- Damaging the power cord may cause fire or electric shock. If the power cord becomes damaged, ask your dealer where the unit was purchased to replace it.
- Unplugging or plugging in the power plug with wet hands may cause electric shock.
- When using this unit, make sure not to cover the AC adapter with a comforter or place objects on top of it. Be sure to observe the following items when using this unit:
 - Do not place any objects around the AC adapter.
 - Do not cover the AC adapter with a cloth.
 - Do not stack two or more AC adapters.
 - Do not put the AC adapter in a box, etc.
- Do not wrap the AC adapter cord or power cord around the AC adapter.

■ On handling

- When carrying the unit, be careful not to damage it or injure yourself.

Important notes

■ Precautions regarding the use of this product

- This product is not intended for treatment, diagnosis, or other medical care.
- Do not use this product for purposes other than odorant release.

■ Precautions regarding odorant injection


- Customers shall be responsible for decision-making and management regarding the handling of the odorants to be injected into the sub-cartridges (sold separately).
- The sub-cartridge material uses anodized aluminum, POM, and PET. Use odorants that do not affect the material.
Furthermore, the material may be subject to change without notice. Customers shall be responsible for making decisions about the compatibility of odorants.
- Sony does not assume any responsibility for the odorants used.
- Follow the steps listed in this Startup Guide or the Help Guide to carry out the series of odorant injection tasks. If the work is not performed correctly, the odorant may leak, adhere to the unit, and leave an odor.

Checking the package contents


In the event that any items are missing, please contact the dealer where the unit was purchased. The numbers in parentheses indicate the number of each item.

The custom cartridge, cartridge clammer, nose guide, and odorant release port cover are optional accessories.


- **Main unit (1)**

- **AC adapter (1)/power cord (2)**

- **Dedicated tablet* set (1)**

The dedicated "Scent Canvas™" app is pre-installed.

- **Startup Guide (this manual)**


* The specifications of the dedicated tablet vary depending on the time of purchase.

EN

Optional accessories


- **Custom cartridge**

(hereinafter, listed as "cartridge")


The cartridge is equipped with 40 odorant injecting sub-cartridges. Cartridge stickers and replacement deodorization filter are included.

- **Cartridge clammer**

Use the clammer together with a commercially available torque wrench to open and close the sub-cartridges. The sub-cartridges can be hermetically sealed.


- **Sub-cartridge set**


This is a set for adding sub-cartridges. Cartridge stickers and replacement deodorization filter are also included.

- **Nose guide/odorant release port cover**

Nose guide


Odorant release port cover

Install a new nose guide on the unit for each user. When not in use, install the odorant release port cover on the unit.

Custom cartridge preparation flow

This manual describes the initial settings in the “Scent Canvas” app, the process of injecting odorant into a sub-cartridge, and the process of linking a custom cartridge (sold separately) with the “Scent Canvas” app, in preparation for using a custom cartridge with the unit.


Preparation Step 1: Arranging all of the necessary items


To inject the odorants into the sub-cartridges and prepare the cartridge, prepare the following items in advance.

Optional accessories

- Custom cartridge

- Cartridge clumper

Commercially available items

The necessary items are not limited to those listed below. Arrange the necessary items according to the odorants to be used and carry out the odorant injection work with the appropriate clothing and environment. If anything is unclear, please contact the dealer where the unit was purchased. The illustrations are for reference only.

- Odorants to be used (maximum of 40 different types)


Note

Sony does not assume any responsibility for the odorants used.


- Torque wrench

Prepare a torque wrench with an insertion angle of 6.3 mm to 6.35 mm (1/4 inch) that can set a torque of 3.5 N·m.


- Micropipette and other dispensing burettes

These are used to inject the odorants. Up to 0.66 mL of odorant can be injected in a sub-cartridge. Use the appropriate product according to the volume and odorant to be injected.

(The illustration shows a micropipette for reference only.)

- Protective glasses


- Protective gloves

- Trash bag for tip disposal


This bag is used for the disposal of used micropipette tips.


Preparation Step 2: Configuring the initial settings with the “Scent Canvas” app

Charge the dedicated tablet to a sufficient level in advance. Check the tablet’s Operating Instructions for the charging method.


- 1 Turn on the dedicated tablet and launch the “Scent Canvas” app.
- 2 Follow the on-screen instructions to select your country/region, check the Software License, and set the language.
When the check finishes, the Home screen is displayed.



- 3 Tap **≡ (Settings)** in the upper right part of the Home screen followed by the **[Admin Mode] (On/Off)** switch to turn Administrator Mode on.

- 4 Tap **≡ (Settings)** in the upper right part of the Home screen and then tap **[Cartridge Settings]**.


- 5 Tap **[+ Add a new setting]** in the upper right part of the screen and then tap **[Create]** on the **[Do you want to create a new cartridge setting?]** screen.
The **[Cartridge Settings]** screen is displayed.

The odorant registration screen is displayed.


- 6 Tap **[Unnamed]** in the upper left part of the screen and enter the cartridge name.

7 Tap [Edit] for the number of the sub-cartridge to register an odorant on the [Cartridge Settings] screen.

8 Enter the odorant information on the [Edit Cartridge] screen.

Enter the [Odorant name] (*Required). And enter the thumbnail image, description, and date of odorant injection as needed.

[Image]:

Save any image to the tablet and set it as the thumbnail image.

[Odorant name] (*Required):

Register the name of the odorant used for each sub-cartridge number (maximum of 24 characters).

[Description]:

Enter the explanatory text for each odorant (maximum of 210 characters).

[Creation date]:

Set the date of odorant injection.

9 Tap [Save] when the odorant information has been entered.


10 Repeat steps 7 through 9 to register the odorants.

Up to 40 different types of odorants can be registered.

11 When all of the odorants have been registered, check the registration details and tap [Complete].

Tap [Preview] to display the cartridge diagram and check the number and position of the sub-cartridges with registered odorants.

Tap [Complete] to save the cartridge settings with the entered details.


Hints

- You can additionally purchase custom cartridges and register multiple custom cartridge settings (maximum of 20).
- You can use a CSV-formatted file to batch register the text for [Odorant name] and [Description] entry. For details about how to use batch registration, refer to the Help Guide.

Preparation Step 3: Removing sub-cartridges from the cartridge

- 1 Place the cartridge on a flat surface.
- 2 Check the sub-cartridge into which the odorant will be injected.



Check the position of the sub-cartridge with the number that registered the odorant with the "Scent Canvas" app.



- 3 Remove the sub-cartridge with the number verified in step 2.

Tilt the sub-cartridge in direction (A) and pull it out.

Sub-cartridges on the outside

Sub-cartridges on the inside


Preparation Step 4: Affixing stickers to the sub-cartridges


1 Write the odorant name and other notes on the stickers included with the cartridge.

Write the sub-cartridge number, odorant name, and other notes while checking the odorant on the Preview screen of the "Scent Canvas" app (refer to step 11 on page 9).

It is recommended that you write with a permanent marker.


2 Affix the sticker in the center of the top side of the sub-cartridge.

Affixing the cartridge name sticker

It is recommended that you write the cartridge name on the rectangular sticker. This makes it easy to identify when using multiple cartridges. Affix the sticker with the cartridge name written on it above the deodorization filter.

Sticker affixing position

Preparation Step 5: Injecting the odorant into the sub-cartridge

Please do the following before starting the odorant injection work.


- Set the torque wrench to 3.5 N·m. For details on how to use the torque wrench, please check the instruction manual for the torque wrench you are using. Opening and closing the sub-cartridge with a torque that is not set to 3.5 N·m may damage the cartridge clamber and sub-cartridge or cause the odorant to leak.
- Set the dispensing volume of the dispensing burette to a maximum of 0.66 mL. For details on how to use the dispensing burette, please check the instruction manual for the dispensing burette you are using.
- Please make preparations suited to the odorant you are using such as wearing protective glasses, protective gloves, etc.

Follow the steps below to carry out the series of odorant injection tasks. If the work is not performed correctly, the odorant may leak, adhere to the unit, and leave an odor.


Opening the sub-cartridge

- 1 Set the sub-cartridge removed from the cartridge in the cartridge clamber (bottom).**


Check that the sub-cartridge does not move once set.

- 2 Set the cartridge clamber (top) over the cartridge clamber (bottom).**


- 3 Use the torque wrench to turn the cartridge clamber (top) counterclockwise and loosen the sub-cartridge (top).**

- 4 Remove the cartridge clamber (top) and take out the sub-cartridge.**

- 5 Manually turn the sub-cartridge (top) counterclockwise to remove it and place the sub-cartridge (bottom) on the table with the top side of the metal cap facing up.**


Injecting the odorant into the sub-cartridge

Use the appropriate dispensing burette according to the volume and odorant to be injected.


Here, the steps are explained using a micropipette for injection.

1 Have the odorant and dispensing burette ready at hand.


2 Attach the tip to the micropipette.

3 Aspirate a maximum of 0.66 mL of the odorant with the micropipette.

4 Inject the odorant from the hole in the center of the top side of the metal cap.


Notes

- Up to 0.66 mL of odorant can be injected. If more than 0.66 mL is injected, the odorant may leak and cause the unit to smell or malfunction.
- In the event that the odorant leaks, wear protective gloves and wipe it off appropriately with a cloth, etc. according to the odorant.


Sealing the sub-cartridge hermetically

1 Place the sub-cartridge (bottom) on a flat surface and manually turn the sub-cartridge (top) clockwise to temporarily tighten it while keeping it level.


2 Keep the sub-cartridge level and set it in the cartridge clamber (bottom).

Check that the sub-cartridge does not move once set.


3 Set the cartridge clamber (top) over the cartridge clamber (bottom).


4 Use the torque wrench set to 3.5 N·m to turn the cartridge clamber (top) clockwise until it clicks and tighten the sub-cartridge (top).

When you hear it click, the sub-cartridge is hermetically sealed. (The sound indicating that it is sealed may differ depending on the torque wrench used.)

5 Remove the sub-cartridge from the cartridge clamber and wait approximately 30 minutes.

It takes approximately 30 minutes for the odorant to permeate inside the metal cap. Do not move the sub-cartridge while you wait.

Place in a level position for approximately 30 minutes.

6 Repeat steps of "Opening the sub-cartridge" (page 12) through "Sealing the sub-cartridge hermetically" (page 13) to inject the odorants into the remaining sub-cartridges.


Notes

- Tightening the sub-cartridge with a torque that is not set to 3.5 N·m may damage the cartridge clamber and sub-cartridge or cause the odorant to leak.
- If the specified torque is not reached while tightening with the torque wrench, the cartridge clamber may be defective. To replace the cartridge clamber, please contact the dealer where the unit was purchased.
- A sub-cartridge can only be used once. Do not reuse a sub-cartridge that was injected with an odorant one time.
- Do not open a sub-cartridge that was injected with an odorant and sealed. Doing so may cause the odorant to leak.

Preparation Step 6: Installing the sub-cartridges in the cartridge


- 1 Place the cartridge on a flat surface.
- 2 Check the install position of the sub-cartridge.

Install the sub-cartridge in the position indicated by the number on the sticker affixed to the sub-cartridge.

- 3 Install the sub-cartridge in the cartridge.



Orient the sub-cartridge so that the vane part is on the upper side and install it directly into the cartridge.

Preparation Step 7: Setting the cartridge into the unit

- 1 Lay the unit down and remove the front cover.


Slide the front cover OPEN lever, turn the front cover slightly in the counterclockwise direction, and remove the front cover.



- 2 Remove the cartridge cover.

Open the 5 buckles fastening the cartridge cover to remove it.


When using the cartridge for the first time, remove the cushioning material from the central part of the cartridge cover.

3 Install the cartridge.


Align the cartridge notch with the protrusion and install the cartridge so that the "SENSOR VALVE" logo appears upright when the underside of the unit is in front of you.

Firmly close the buckles. If the front cover is attached without firmly closing the buckles, the cartridge cover may be damaged.


5 Attach the front cover.

Turn the front cover clockwise until it clicks to attach it.

4 Attach the cartridge cover.


With the "SENSOR VALVE" logo appearing upright when the underside of the unit is in front of you, attach the cartridge cover so that the protrusion passes through the opening and fasten the 5 buckles.


Preparation Step 8: Turning on the power to the unit

1 Connect the power cord of the AC adapter to the unit and the outlet.

When connecting the power cord to the AC adapter, select the power cord with the plug that matches the shape of the outlet.


2 Turn on the power switch on the right side.


Preparation Step 9: Connecting the unit and tablet

Charge the dedicated tablet to a sufficient level in advance. Check the tablet's Operating Instructions for the charging method.

- 1 Turn on the dedicated tablet and launch the "Scent Canvas" app.
- 2 Tap  (Connection) in the upper left part of the app screen.

3 Tap [Select a device].

4 Choose and tap the unit to connect from the displayed list.

The unit and dedicated tablet are connected via Bluetooth® communication.


5 When [Connected] is displayed, tap X (Close) in the upper right part of the screen.

Preparation Step 10: Linking the cartridge and the cartridge setting


Before linking the cartridge, check that the sub-cartridges with the injected odorants are correctly installed in the cartridge. For details, see "Preparation Step 6: Installing the sub-cartridges in the cartridge" (page 15).

1 Tap [+ Link to the cartridge] on the [Cartridge Settings] screen.


2 Check that the linking is complete.

When the linking is complete, the cartridge ID is displayed in the status display area in the upper left part of the screen, and the indicator changes from ! (Cartridge unconnected) to ✓ (Cartridge connected).

3 Tap [←] (Return).

■ Linking with other cartridge settings


If multiple cartridge settings are registered, you can choose which cartridge setting to link to the cartridge installed in the unit. For details on linking with other cartridge settings, refer to the Help Guide.

Overview of the four types of templates

The “Scent Canvas” app provides the following four types of templates.


You can use the templates to customize the setting items or to add a description and create a process tailored to the application (setting file that summarizes the odorant release conditions).

For more details about how to use the templates, refer to the Help Guide.

- **Index**
Select any odorant from a list of registered odorants to release an odorant.
- **Random**
This type of template randomly releases an odorant before offering choices and having the user answer. Explanations for each question and the number of correct answers for the question can be displayed.
- **Preset (single)**
This type of template releases one arbitrary odorant before offering choices and having the user answer. Explanations for each question and the number of correct answers for the question can be displayed.
- **Preset (multiple)**
This type of template continuously releases multiple odorants before offering choices and having the user answer. Explanations for each question and the number of correct answers for the question can be displayed.

Help Guide (Web manual)

The Help Guide is an online manual that can be viewed via the Internet using a computer or smartphone.

For information about how to create a process using the template or how to use the unit with the custom cartridge, refer to the Help Guide.

https://rd1.sony.net/help/nos/nosdx1000/h_zz/

Troubleshooting

If the unit or “Scent Canvas” app does not work properly, follow the steps below.

- ① **Save the data recorded by the “Scent Canvas” app.**
Exit the process with the “Scent Canvas” app.
- ② **Turn the power to the unit off and on.**
Press the power switch on the right side of the unit to turn off the power and then turn it on again.
- ③ **Relaunch the “Scent Canvas” app.**
 1. Press and hold the “Scent Canvas” app icon on the tablet Home screen and tap [Ⓢ] – [Force stop] – [OK]. (Depending on your tablet, the steps may be: press and hold the “Scent Canvas” app icon, and select [App info] – [FORCE STOP] – [OK].)
 2. Return to the Home screen and launch the “Scent Canvas” app.

Check the support website.

<https://www.sony.net/nos-support-global/>


Contact the dealer where the unit was purchased.

Precautions

■ Use of this product

- Avoid using this product in a location with strong odors.
- If the released odorants remain inside the room, ventilate the room.
- Be sure to use authorized products for custom cartridges, sub-cartridges, cartridge clampers, nose guides, and odorant release port covers, which are sold separately.
- Do not use any adapters other than the supplied AC adapter.
- Do not block the exhaust port on the rear of the unit. If the exhaust port is blocked, it will become unable to deodorize the unit and may stop working correctly.

■ Handling of personal information

- Personal information, answer results, and other data recorded by the "Scent Canvas" app are stored on the dedicated tablet. Please handle such data with sufficient care.
- The personal information, answer results, and other data recorded by the "Scent Canvas" app are not encrypted. Please beware of data leakage via the USB port of the dedicated tablet or unauthorized access by third parties over the Internet.

■ Heat generation

The dedicated tablet may heat up during operation. Wait until the temperature decreases before using the tablet.


■ Handling custom cartridges (sold separately) and sub-cartridge sets (sold separately)

- When injecting odorant into a sub-cartridge, perform the work by following the procedures listed in the Help Guide or the Startup Guide (this manual).
- Up to 0.66 mL of odorant can be injected. If more than 0.66 mL is injected, the odorant may leak and cause the unit to smell or malfunction.
- In the event that the odorant leaks, wear protective gloves and wipe it off appropriately with a cloth, etc. according to the odorant.

- A sub-cartridge can only be used once. Do not reuse a sub-cartridge that was injected with an odorant one time.
- Do not subject the cartridges or sub-cartridges to strong impact such as dropping them on the floor. Doing so may cause them to malfunction.

■ Maintenance

- Carefully wipe away any dirt or moisture on the unit before placing it in storage. If the unit is very dirty, wipe it clean with dehydrated ethanol.
- If you detect other odorants, clean the front cover odorant release port and the cartridge cover odorant outlet (both the front and back sides) with a cotton swab moistened with dehydrated ethanol. If you are concerned about the odorant, repeatedly clean these areas with dehydrated ethanol several times. Depending on the odorants used, take appropriate measures such as wearing protective gloves, etc.

■ Storage

- To prevent contamination and deterioration, do not store the unit or nose guides, odorant release port covers, custom cartridges, sub-cartridges, or cartridge clampers, which are separately purchased accessories, in locations with high temperatures, direct sunlight, high humidity, strong odors, etc.
- Store cartridges equipped with sub-cartridges that have injected odorants under conditions that are suitable for the odorants used.

■ Notes on using wireless LAN

Please note that Sony cannot make any guarantees regarding malfunction, failure, or damage to the dedicated tablet caused by connecting to the Internet, except when instructed by the manufacturer or dealer.

Disposal

■ Disposal of the unit, nose guide, odorant release port cover, and cartridge clamper

Discard these items according to the ordinances, etc. of the region or municipality where you live.

■ Disposal of the cartridge and sub-cartridges

Discard these items according to the ordinances, etc. of the region or municipality where you live.

If you have any problems concerning disposal, please contact the dealer where the unit was purchased.

■ Disposal of the dedicated tablet

Follow the method listed in the tablet's Operating Instructions to dispose of the dedicated tablet.

Specifications

Main unit

Model No.

NOS-DX1000

Maximum external dimensions (W × H × D)

Approx. 37.4 cm × 37.8 cm × 25.3 cm

Mass

Approx. 5.2 kg

Power source

DC 12 V (using the supplied AC adapter)

Power consumption (maximum)

Approx. 6.0 W

Eco and energy saving

Power consumption (Off mode): 0.3 W

Power consumption (Standby mode/Low power mode): Not available

Power consumption (Networked standby mode): 2.0 W (Bluetooth)

Networked standby mode activated: Within 20 minutes (default)

Recommended environment

Ambient temperature: 20°C to 25°C

Ambient humidity: 40% to 70%

Odorant source

Cartridge system: a maximum of 40 types of odorants can be loaded into one cartridge

Odorant switchover time

Approx. 3 seconds

Odorant release time

1 to 10 seconds

Odorant release method

Gaseous diffusion method (non-heated, dry air method)

Operating method

Operate from the "Scent Canvas" app on the supplied dedicated tablet

Communication method

Bluetooth Ver. 5.0 (Bluetooth Low Energy)

Dedicated tablet

Please see the Operating Instructions for the tablet.

Trademarks and copyrights

- Tensor Valve is a trademark or registered trademark of Sony Group Corporation.
- Scent Canvas is a trademark or registered trademark of Sony Group Corporation.
- The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Sony Group Corporation and its subsidiaries is under license.
- Wi-Fi®, Wi-Fi Alliance®, and Wi-Fi Direct® are registered trademarks of Wi-Fi Alliance.
- Android is a trademark or registered trademark of Google Inc.

- Please note that the specifications and external appearance may be subject to change without notice for improvement.
- The actual values of specifications may deviate from the estimated ones depending on the actual environment and method of use.

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- Tensor Valve è un marchio o marchio registrato di Sony Group Corporation.
- Scent Canvas è un marchio o marchio registrato di Sony Group Corporation.
- Android è un marchio o marchio registrato di Google Inc.